



Performance Meets Compliance



AMSCO® 400 Series
Medium Steam Sterilizer



Performance meets compliance:

The AMSCO® 400 Series Medium Steam Sterilizer is designed for performance in today's sterile processing environments and is AAMI ST8:2008 compliant.

Built to process 47% more instruments by weight*; the AMSCO 400 Series Medium Steam Sterilizer is validated to safely sterilize heat and moisture stable materials in packs up to 25 pounds. This allows you to efficiently handle high volume sterilization and heavier trays, while meeting today's stricter compliance guidelines.

- AAMI ST8:2008 compliant
- Processes up to 25 pound packs and trays
- Keeps same footprint as previous models



**than previous generation*

Effortless operation

Easy to open and easy to load, the AMSCO 400 Series Medium Steam Sterilizer will help maintain a safe and productive work environment. Whether you prefer a power or manual door, opening requires minimal effort.

- Easy to open doors
 - Hinged doors feature a patented single hand, low-effort handle locking mechanism
 - Power sliding doors are easily controlled with a touch of a button
- Loading car and transfer carriage
 - Unique docking station simplifies loading and unloading



Quick, efficient control

Equipped with an intuitive, color, touch-screen control panel, the AMSCO 400 Series Medium Steam Sterilizer allows your staff to be more efficient with simple cycle selection.

- Comes standard with pre-programmed cycles including prevac, gravity, liquid, leak test and DART – all validated to AAMI ST8:2008 standards
- Up to 12 programmable cycles
- One-touch operation
- Customizable cycle options
- Eye-level, ergonomic control
- Control countdown screen visible from 15'



Convenient one-touch control at eye level.



Cycle countdown visible from 15'

Safe, strong and green

- Non-lubricated, steam-activated door gasket ensures no air leaks
- Door stops automatically if there is an obstruction
- Fifteen-year chamber warranty
- Elliptical or round chamber, jacket and door constructed with 316L stainless steel
- Automatic utilities start-up and shut-down
- Electronic water-saving device



AMSCO® 400 Series Medium Sterilizer with hinged manual door
26" x 37.5" (660 x 950 mm)



AMSCO® 400 Series Medium Sterilizer with vertical power door
26" x 26" (660 x 660 mm)



AMSCO® 400 Series Medium Sterilizer with horizontal power door
26" x 37.5" (660 x 950 mm)

STERIS offers additional solutions that fit perfectly into your SPD. Promote productivity, traceability and compliance with our workflow management system and ProConnect® Remote Monitoring. Talk to your STERIS representative to learn more.

